


# i400 – Clinical excellence meets design-led style



The i400 combines the most advanced iPulse™ technology with an exclusive design for premium skincare establishments. This latest addition to the celebrated iPulse™ range perfectly combines aesthetics and clinical excellence.

# i400



## Treat large areas quicker

Like all iPulse™ systems, the i400 is designed to produce the most effective clinical results in the quickest time, at the lowest possible cost to the practitioner. The system's large 8.9 cm<sup>2</sup> beam of light means areas of the body are treated quicker with the i400, and the perfect level of energy generated by icePulse, a proprietary power supply, makes each pulse of light more effective, resulting in fewer treatments. This saves the practitioner time and money.

## Lighter and sleeker

Cosmetic devices using light energy have traditionally been heavy, bulky machines, while smaller devices (sometimes known as mini-IPLs) compromise on clinical results. The advanced icePulse technology behind the i400 allows it to deliver the highest clinical performance from a small 31.8 lb. unit. This compact unit is portable within the treatment environment and fits neatly onto a mobile cart or standard desktop.

## Two options

The i400 is available in two versions. One is a small table-top unit and the other comes with an integral cart designed to suit a clinical environment where the device needs to be used in different areas.



# i400 – an IPL for the style conscious




## Technical Specification

Light source:	twin Xenon lamps
Wavelength range:	530 – 1200 nm
Spot size on tissue:	8.9 cm <sup>2</sup>
Pulse duration:	10 – 110 ms
Pulse sequence:	single or triple pulse flash
Triple pulse:	5 – 20 ms on-time 7 – 30 ms off-time
Repetition rate:	average of 0.25 Hz
Energy density (fluence) range:	0 – 20 J/cm <sup>2</sup>
Cooling:	fan, air-cooling
Dimensions for table-top:	7.1" (high) 19.7" (wide) 15.6" (deep)
Dimensions for cart version:	37.8" (high) 19.7" (wide) 15.6" (deep)
Weight:	31.8 lbs. (table-top) 43.0 lbs. (with cart)

## Improved User-Interface

Adjustment of settings has been simplified so the device can be easily set up for different clients. The practitioner selects the treatment parameters based on a combination of skin types and desired results to provide exactly the right treatment for each individual client.



Touch sensitive control and LCD display allow for easy set-up.

## Ready to use

The i400 plugs into a standard wall outlet and is ready to use immediately. Lamp cartridges provide the optimum light energy for treatment throughout their fixed lifetime; they are just as powerful on their 10,000<sup>th</sup> shot as their first. The computer alerts the operator when the lamp needs to be replaced... a ten second task for the operator just like changing a normal light bulb.



The i400 is available in two versions—table-top and standalone with integral cart.

## All Inclusive Package

The i400 system features:

- One multi-procedural hand piece
- One lamp assembly
- Two pairs of safety glasses
- One user's manual
- One year parts and labor warranty
- FDA cleared for multiple indications



San Diego | Boston | London | Sydney

**Med-Aesthetic Solutions, Inc.**  
2033 San Elijo Ave. #200  
Cardiff by the Sea, CA 92007  
Customer Support: 877.733.7627  
505.341.2577 | Fax: 505.341.2575  
[www.medaestheticsolutions.com](http://www.medaestheticsolutions.com)

iPulse is a registered trademark and patented technology of CyDen Limited © 2008, Med-Aesthetic Solutions, Inc. All rights reserved. Med-Aesthetic Solutions, Inc. and its logo are trademarks of Med-Aesthetic Solutions, Inc.

